

Check-list

- Check connections
 - Laptop has power
 - Kinect has power
 - Kinect is plugged into a fast USB port
 - Kinect Studio is running
 - Kinect Studio is connected to the Kinect
 - Kinect's infrared lights are visibly on
 - Leap is plugged into a fast USB port
 - Leap's infrared lights are visibly on
 - Myo USB is plugged in
 - Myo is orientated with the battery light towards the wrist
 - Myo is calibrated and sync'd for participant
- Run session
- Move participant data into folder for preprocessing